

Specialist Wound Care Practitioner Dressing Formulary

(To be requested by a clinician (Oxleas District Nurses) only after discussion with CWCT Nurse Specialist)

Dressing	Category	Sizes	Pack size
Acticoat	Primary antimicrobial agent for infected wounds	5 x 5cm 10 x 10cm	5 5
Acticoat flex 3	Primary antimicrobial agent for infected wounds – can be used under TNP	10 x 10cm	5
Actilymph hosiery	Diagnosed lymphoedema with venous component	Small, Medium Large, X-large XX-large	1pair Rx 2 pairs up to 6 months
Anabact	Primary wound gel containing Metronidazole for topical application – for malodorous, infected, necrotic wounds	15g 30g	1 tube 1 tube
Maggots	Larvae to debride recalcitrant slough Biobag or Free-range	Various sizes (Unlicensed, not listed in DT)	NA
Mepitel	Soft silicone contact layer	5 x 7cm 8 x 10cm 12 x 15cm	5 5 5
PICO - Negative pressure dressing kits	Provides 7 days continuous negative pressure step - down from TNP or for cavity wounds low exudate Pump + 2 dressings	10 x 20cm 10 x 30cm 15 x 15cm 15 x 20cm	1 1 1 1
Kerraped Plus boots	Patients with neuropathic foot ulcers or who cannot wear their own footwear because of bulky bandages. Can also be used to offload pressure from the front of the foot.	Small, Medium, Large and Extra large	1pair
Urgotul Silver	Used primarily for plastic wounds and burns but also where silver is needed and aquacel AG not appropriate.	10x12cm	1

NOTE:

- Bank Nurses and Agency Nurses must **NOT** be involved in ordering from Specialist Formulary.
- Order for a maximum 2 weeks .If the wound fails to improve the patient needs to be referred to the CWCT for reassessment.